

**On behalf of Vision Expo, we sincerely thank you for being with us this year.**

**Vision Expo Has Gone Green!**

We've eliminated paper session evaluation forms. Please be sure to complete your electronic session evaluations online when you login to request your CE Letter for each course you attended!

Your feedback is important to us as our Education Planning Committee considers content and speakers for future meetings to provide you with the best education possible.



1

---

---

---

---

---

---

---

## Financial Disclosure Statement

**Andrew Bruce has received honorarium from:**

- VSP Optics/UUniversity
- Mitsui Chemicals America, Inc.
- Optical Training Institute

- All relevant relationships have been mitigated
- He has NO financial interest in any product presented in this course.

2

---

---

---

---

---

---

---

## Index Matching For Optical Excellence

ABO Approved Level II

Course Written by Andrew S. Bruce  
LDO, ABOM, NCLEI, FCLSA

3

---

---

---

---

---

---

---

## Objectives

Upon completion, the participant should be able to...

Discuss refractive index and its significance, with regards to both lens materials and coatings.

Communicate the benefits that result when the refractive indices of the hard coating and lens substrate are closely matched.

Appreciate the benefits of recommending MR™ premium high index.

4

4

---

---

---

---

---

---

---



5

---

---

---

---


---

---

---

## Which Material/Substrate Is Best?

Depends on...



6

6

---

---


---

---

---

---

---



## Refractive Index

Refractive Index,  $n = \frac{\text{Speed of light in air, } c}{\text{Speed of light in second medium, } v}$

**Refractive Index:** Ratio between speed of light in air to speed of light in a second medium

- A measure of a medium's refractive properties
- Higher refractive indices = thinner lens
- Is all high index created equal?

7

---

---

---


---

---

---

---


---



## Not All High Index is Created Equal

Mitsui Chemicals is the world's major producer of the chemistry used to make the MR™ series of thiourethane-based, premium high index

- MR™ high index includes:
  - **MR-8™** (1.60 / Abbe 41)
  - **MR-7™ & MR-10™** (1.67)
  - **MR-174™** (1.74).



8

---

---

---

---

---

---

---

---



Technologies

[sdctech.com](http://sdctech.com)

38 YEARS OF PRODUCT INNOVATION

**CrystalCoat™**



9

---

---

---

---


---

---

---

---

### Benefits of MR™ Thiourethane Premium High-Index



- Thin, light-weight lenses
- High Abbe values / heat distortion temp.
- Easy to tint / 100% UV attenuation
- Excellent mechanical strength
- Little to no discoloration, over time
- Superior coating adhesion/durability
- Index matched hard coatings

10

---

---

---

---

---


---

---

---

### Acrylic High Index

- Many incorporate hard coatings index-matched for low index resins, **NOT** high index substrates
- Creates an index mis-match
- Light perceives the coating-substrate interface as two separate entities
- Results in less-than-ideal outcomes.



11

---

---

---

---

---

---

---

---

### Key Benefits of Index Matching

#### 1: Minimized Reflectance

- Light perceives the coating-substrate interface as a single entity
- The closer the index matching, the lower reflectance

#### 2: Strong Coating Adhesion

- Demonstrates superior durability in industry standard testing to withstand daily wear and tear

#### 3: Minimal Birefringence

- Perceived single interface results in minimal birefringence
- Delivers optimal clarity and cosmetics.

12

---

---

---

---


---

---

---

---

## 1: Minimized Reflectance



- A perfect index match would result in a flat line or zero reflectance
- Perfect index match not practically possible
- The closer the match, the lower reflectance, hence, a smoother line.

---

---

---

---

---

---


---

---

13

## 2: Coating Adhesion Testing with MR-8™

### Primary Adhesion Test



14

---

---

---

---

---


---

---

---

## 2: Coating Adhesion Testing with MR-8™

### Boiling Water Test



15

---

---

---

---

---

---

---

---

15

## 2: Coating Performance Test Results

Hard Coating Refractive Index	1.60	1.66
Primary Adhesion Test	100%	100%
Boling Water Test (After 30 Minutes)	100%	100%
Bayer Test	3.7 (on MR-8™)	2.7 (on MR-7™)



- Bayer Test & Bayer Ratio:
- Measures haze gain on coated lens vs. uncoated CR39 reference lens
  - Bayer ratio is a measure of the scratch resistance of a coated lens vs. uncoated reference lens.

16

---

---

---

---

---


---

---

---

## 3: Minimal Birefringence

### Mis-Matched Hard Coating



Significant birefringence

### Index-Matched Hard Coating



Minimal birefringence

17

---

---

---

---

---

---

---

---

How Do We  
Index Match?

18

---

---

---

---

---


---

---

---

# How Do We “Index Match”?

Standard 1.49 index  
hard coating is  
silica based



19

---

---

---

---

---

---


---

---

# How Do We “Index Match”?

Standard 1.49 index  
hard coating is  
silica based

To increase its index,  
metal oxides are  
added to coating resin



20

---

---

---

---

---

---

---


---

# How Do We “Index Match”?

Standard 1.49 index  
hard coating is  
silica based

To increase its index,  
metal oxides are  
added to coating resin

Mostly used oxides:  
Titania dioxide or  
Zirconia dioxide



21

---

---

---

---

---

---

---

---

7


## How Do We “Index Match”?

Standard 1.49 index  
hard coating is  
silica based

To increase its index,  
metal oxides are  
added to coating resin

Mostly used oxides:  
Titania dioxide or  
Zirconia dioxide

Precisely controlled  
addition prevents  
*nanoparticle*  
*agglomeration*



22

---

---

---

---

---


---

---

---

## Nanoparticle Agglomeration

- Clumping together of metal oxide nanoparticles
- Can result in degraded quality and “spotty” lenses
- MRT™ series hard coatings comply with strict QC standards
- Ensure great performing lenses.



23

---

---

---

---

---

---

---

---

**CrystalCoat™**



24

## Application Requirements

MRT™ substrates require:

- Contamination-free environment
- Primer layer or caustic etching under optimized conditions
- Sufficient curing to ensure complete caustic removal from the lens.

---

---

---

---


---

---

---

---





25

---

---

---

---

---

---

---

---

### Primary Roles Of A Hard Coating?

- Abrasion resistance
- Stable substrate for AR coating.

26

---

---

---

---

---

---

---

---

### Why Recommend AR Coated Lenses?

27

---

---

---

---

---

---

---

---

## Benefits of Premium AR Coatings

Multiple layers  
minimize lens  
surface reflections

Optimize light  
transmission/  
visual acuity

Enhanced  
durability/long-  
term performance

Complement  
digital lens benefits

Easy maintenance  
and excellent  
warranties

Increased patient  
satisfaction

28

---

---

---

---

---


---

---

---

28

## How AR Coatings “Do Their Thing”



AR layer thickness MUST be precisely  $\frac{1}{4}$  wavelength of the incident light

29

---

---

---

---


---

---

---

---

## What About Premium AR Coatings



30

---

---

---

---

---

---

---

---



31

---

---

---

---

---

---

---



32

---

---


---

---

---

---

---



33

---

---

---

---

---


---

---

34

## Pop Quiz Time . . .

- During the primary adhesion test, a razor blade is used to create cuts through a lens surface at 90° to each other. What do we call this cut area?
- How is a lens material's refractive index increased?
- Who is the world's major producer of the chemistry used to make MR premium high index lens materials?



34

---

---

---


---

---

---

---

---



## Any Questions?

35

---

---

---

---

---

---

---


---

**On behalf of Vision Expo, we sincerely thank you for being with us this year.**

**Vision Expo Has Gone Green!**

We have eliminated all paper session evaluation forms. Please be sure to complete your electronic session evaluations online when you login to request your CE Letter for each course you attended!

Your feedback is important to us as our Education Planning Committee considers content and speakers for future meetings to provide you with the best education possible.



36

---

---

---

---

---

---

---

---

12



### Thank you

For more information, visit the  
MR™ High-Index Materials  
LinkedIn Page:



This ABO Approved Course is supported by an  
educational grant from Mitsui Chemicals.

37

---

---

---

---

---

---

---