


# Compensated Powers (and other ophthalmic conundrums)

**Pete Hanlin, ABOM**  
Sr. Director Professional Solutions  
Essilor of America




---

---

---

---


---

---

---

---

### Same Lens- Different Results




---

---

---

---


---

---

---

---

### PALs are not measured @ OC




---

---

---


---

---

---

---

---




---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

▶ Measured vs. Worn Conditions

ANGLE (°)	SPHERE (D.)	CYLINDER (D.)	AXIS (°)
0	4.00	0.00	-
3	3.99	0.01	90
6	3.94	0.04	90
9	3.87	0.10	90
12	3.77	0.17	90
15	3.64	0.27	90

Measured focimeter powers of an as-worn +4D. SV calculated lens taking into account a tilt of ANGLE.

Has the power really been "compensated" (changed)?  
 No. The lens is still a +4.00 sph.

"Compensated" power indicates how the lensometer will see the lens.

---

---

---

---

---

---

---

---

### Methods of Measurement

**+5.00 -4.00 x 090 FOA vs IOA (CV)**


	Min. sphere (D.)	Max. sphere (D.)
FOA	1.03	5.06
IOA	0.99	4.98
Discrepancies (FOA - IOA)	<b>+0.04</b>	<b>+0.08</b>

**ADD measurement discrepancy FOA (CV v CX)**

Mean sphere (D.)	Far vision	Convex	Concave
		3.05	2.92
	Near vision	5.72	5.40
	Addition	2.48	2.67
Discrepancies (Concave - Convex)		<b>+0.19</b>	

**ADD measurement discrepancy IOA (CX v Wearer)**

Mean sphere (D.)	Far vision	Convex	Wearer
		2.95	3.05
	Near vision	5.46	5.90
	Addition	2.51	2.84
Discrepancies (Wearer - Convex)		<b>+0.33</b>	




---

---

---

---

---

---


---

---

---

---

### Measured power ≠ ordered power



Wearer prescription →

Expected focimeter measurements →

---

---

---

---

---

---

---

---

---

---

### Compensated Powers

**“Compensation” does NOT change the ordered power...**

**...compensation indicates how the ordered power will be seen by the lensometer!**

---

---

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

### Abbe Value and Chromatic Aberration

**Abbe Facts (Fabbies ☺)**

- Plano lenses produce ZERO chromatic aberration
- No LCA at the optical center- regardless of power (if blur occurs looking straight ahead, its not LCA)
- Human eye has >1.00 diopter of ACA
- Studies show >0.12D CA may be noticeable to some (>3.75D of power/prism is required for 0.12D CA)
- If patient non-adapts to poly, but accepts 1.67... the problem was NOT chromatic aberration
- AR does NOTHING to reduce chromatic aberration

---

---

---

---

---

---

---

---

### Abbe Value and Chromatic Aberration

**Abbe Facts (Fabbies ☺)**

- Studies show >0.12D CA may be noticeable to some

**CR-39 (58) = 0.06Δ LCA = 20/21**  
**Trivex (44) = 0.09Δ LCA = 20/22**  
**Polycarb (30) = 0.13Δ LCA = 20/23**  
**1.60 MR-8 (41) = 0.09Δ LCA = 20/22**  
**1.67 MR-7 (32) = 0.12Δ LCA = 20/23**

$LCA = \frac{3.75}{ABBE} \times PRISM$

**Chromatic Aberration is NOT an issue in ANY material (for 93.7% of wearers at least)!**

LATERAL CA	VISUAL ACUITY
0.05 Δ	20/21
0.10 Δ	20/22
0.15 Δ	20/24
0.20 Δ	20/26
0.25 Δ	20/28
0.30 Δ	20/31
0.35 Δ	20/34
0.40 Δ	20/39
0.45 Δ	20/44
0.50 Δ	20/51
0.55 Δ	20/60
0.60 Δ	20/75

Effect of Chromatic Dispersion of a Lens on Visual Acuity. Meislin, D. & Obrecht, G. Am. J. of Optom. & Physiol. Optic. 65:25-26, 1988.

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

---

---

### 5.1 General

Both uncut and edged finished lenses shall meet the following requirements. For lenses produced with compensations to account for as worn correction, the tolerances in the tables in clause 5 apply to those values specified by the manufacturer and not to the prescribed RX.

#### 5.1.1 Distance Refractive Power (Back Vertex Power)

##### 5.1.1.1 Single Vision and Multifocal Lenses

Table 1 – Tolerance on Distance Refractive Power (Single-Vision and Multifocal Lenses)

Sphere Meridian Power	Tolerance on Sphere Meridian Power	Cylinder ≥ 0.00 D ≤ -2.00 D	Cylinder > -2.00 D ≤ -4.50 D	Cylinder > -4.50 D
From -6.50 D to +6.50 D	± 0.13 D	± 0.13 D	± 0.15 D	± 4%
Stronger than ± 6.50 D	± 2%	± 0.13 D	± 0.15 D	± 4%

##### 5.1.1.2 Progressive Addition Lenses

Table 2 – Tolerance on Distance Refractive Power (Progressive Addition Lenses)

Sphere Meridian Power	Tolerance on Sphere Meridian Power	Cylinder ≥ 0.00 D ≤ -2.00 D	Cylinder > -2.00 D ≤ -3.50 D	Cylinder > -3.50 D
From -8.00 D to +8.00 D	± 0.16 D	± 0.16 D	± 0.18 D	± 5%
Stronger than ± 8.00 D	± 2%	± 0.16 D	± 0.18 D	± 5%

---

---

---

---

---

---

---

---

---

---




---

---

---

---

---

---

---

---

---

---